- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению
- zh** 使用说明书

KaWe MedCharge® 4000

QM-1-170C

de	Ladestation KaWe MedCharge® 4000	4
en	Charging Station KaWe MedCharge® 4000	8
fr	Chargeur KaWe MedCharge® 4000	12
it	Caricatore KaWe MedCharge® 4000	16
es	Cargador KaWe MedCharge® 4000	20
pt	Estação de carga KaWe MedCharge® 4000	24
ru	Зарядное устройство KaWe MedCharge® 4000	28
zh	充电器 KaWe MedCharge® 4000	32

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User's Manual Charging Station KaWe MedCharge® 4000

Dear customer, thank you for choosing a KaWe product. Our products are known for their high quality and durability. This KaWe product complies with all EMV and low voltage guidelines.

Please read these instructions thoroughly and carefully before attempting to use this product and heed the maintenance instructions.

Familiarise yourself thoroughly with the operation of this device before using it.

Utilisation: The MedCharge® 4000 Li-Ion charging station is only to be used by authorised and trained personnel. The MedCharge® 4000 Li-Ion charging station charges only 2.5 V KaWe rechargeable battery handles and/or 3.5 V (NiMH) rechargeable batteries as well as 3.5 V Li-Ion rechargeable batteries. The battery type is detected automatically, i.e. NiMH or Li-Ion and the charging characteristics are controlled as is appropriate for each type. Charging cell type Sub C Li-Ion 3.5 V with 2250 mAh; type AA NiMH 2.5 V with 1600 mAh or type AA NiMH 3.5 V with 700 mAh.

Purpose: The rechargeable batteries are to be used as an energy source for the relevant KaWe products.

Unsuited use/contraindication of the products: Any use other than that described here is not in accordance with the intended use of this product. The manufacturer is not liable for any resulting damages. The user alone bears the risk.

Further instructions, maintenance, storage: When used and stored properly, your KaWe MedCharge® 4000 charging station will serve you reliably for many years.

Safety instructions: **1.** Before the first use, the rechargeable handles must be charged at least 10 hours! Caution: do not use dry cell batteries – risk of explosion! **2.** Before

plugging the device in, ensure that the voltage from the outlet is within the permissible range (see technical data) **3.** After charging has been completed, please remove the batteries from the device and unplug the charging station's power cord from the wall outlet. **4.** Only use in closed and dry rooms.

Guarantee: When used under normal circumstances and with attention to these instructions (the User's Manual) we guarantee this product for two years after the date of purchase. Should you need further information or should your instrument require repair, please contact your dealer.

Operation: After the transformer plug has been plugged into the outlet, the charging station is ready for use. (The green LED display will read: "stand by".) When the rechargeable handles are inserted into the station, the LEDs on the respective charging slots will display the charging progress (red = charging, green = fully charged, blinking red = error). After the target charging voltage has been reached, the LED changes from red to green, i.e. the rechargeable handle has been fully charged. Charging begins automatically as soon as the handle is inserted into the charging slot. The charging slot can be adapted as needed to accommodate the different types of charging handles. The rechargeable battery can be removed from the handle / replaced by unscrewing the cap over the battery compartment.

Performance characteristics:

- Input voltage 100V-240 V +/- 10%, 50-60 Hz
- Max. power 1.2 A
- Efficiency factor 90-93%
- Power factor >0.99
- Charging cells NiMH und Li-Ion
- Nr. of cells NiMH 2-3 cells, Li-Ion 1 cell
- Chargeable capacity min. 250 mAh, max. 4500 mAh
- Input voltage charging station 9 VDC +/- 10%
- Power consumption in standby < 0.5 W
- Rated power consumption approx. 4.5 W
- Output voltage in standby < 9VDC

- Rated output voltage2 x 3.6 VDC
- Rated charging current.....2 x 250 mA +/- 10%
- Charge end voltage.....6 VDC +/- 1%
- Charge end power10-30 mA for Li-Ion
- Return current without power< 0.1 mA at 4.2 V
- Sustained charging powerNiMH 3-8 mA,
Li-Ion has no sustained charging power
- Electric strength4 kV
- Safety ratingII
- Operating temperature0°C - +40°C
- Storage temperature.....-25°C - +70°C

Tests: CE, Innova GS, cCSAus, SAA

Norms: IEC/EN 60335-1, IEC/EN 60335-2-29, UL1310

Safety features: **1.** Protected electronically against short circuiting, protected mechanically against reverse polarity. **2.** Protected against overheating (linear derating). **3.** Built-in reverse voltage protection (safety). **4.** Automatic shut-off with grid overvoltage. **5.** Automatic temperature monitoring of the rechargeable batteries.

Caution: **1.** Please use only accessories recommended by KaWe. Any other connection could cause a fire, electrocution or injury. **2.** The AC transformer that comes with the product produces 9 V direct voltage and 1000 mA direct current. Please do not use any other transformer that does not meet these requirements as this can damage the charger or the transformer. **3.** Please do not alter the plug that comes with the product. Should the plug not be compatible with the outlet, have the proper outlet installed by a qualified electrician. Improper connection could cause electrocutions. **4.** Do not use the charging unit if it has been damaged. Please have it repaired by a qualified technician. **5.** Please do not disassemble the charger. Incorrect reassembly could cause a fire or electrocution. **6.** During the charging process, the batteries and the charger can become warm. **7.** Please never use the charger as a source of electricity for other electrical appliances. **8.** Protect the charger and the transformer from moisture.

Cleaning: The outside surface of the charging station can be cleaned with a damp, soft, lint free cloth. To disinfect the unit, dampen the cloth with alcohol. The unit may not be immersed in any solution or cleaning agent!

Optional accessories: Mounting bracket for attaching the unit to the wall.

	At 30% - 75% rel. humidity
Storage	+5°C to +45°C
Transport	-20°C to +50°C
Operation	+5°C to +40°C

Disposal: Please dispose of your electronic waste in accordance with legal regulations. Do not dispose of these units in regular household garbage, but at a recycling centre that accepts electronics.

Common accessories: Further information about these products can be obtained in our Internet shop at: www.kawemed.de.

Manufacturer: **KaWe**

Contacting the manufacturer: Address or telephone number of the dealer or dial +49-7141-68188-0.

Symbol key:

	Manufacturer		Separate disposal of electric and electronic devices
	Date of manufacture		Disposal of batteries
	Heed the User's Manual		Complies with relevant EU guidelines
	Caution		

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KIRCHNER & WILHELM GmbH + Co. KG